

GENEMER SICKLE CELL

Background

Sickle cell anemia is an autosomal recessive disease. A single base change (A to T) in the β -globin chain causes the substitution of amino acid glutamine to valine: the cause of the disorder sickle cell anemia. The resulting mutant globin chain is termed as the Hb S. Hemoglobin S is freely soluble when fully oxygenated, under conditions of low oxygen tension the red cells become grossly abnormal assuming a sickle shape leading to aggregation and hemolysis. Homozygous Hb S is a serious hemoglobinopathy found almost exclusively in the black population. About 8% of American Blacks are carriers and about 0.2% are affected. Heterozygotes (sickle cell trait) are clinically normal, although their red cells will sickle when subjected to very low oxygen pressure in vitro.

DNA analysis for the sickle cell mutation is done by specific amplification of the DNA region spanning the mutation using polymerase chain reaction followed by enzymatic cleavage of the amplified product. Sickle cell mutation abolishes a restriction endonuclease site (Dde I). Electrophoretic resolution of the fragment pattern reveals the presence or absence of the mutation. Clear genotyping of normal, carrier and homozygous DNA is achieved.

Protocol for Sickle Cell DNA Genotyping

DNA EXTRACTION FROM BLOOD

1. Collect 0.5mls of whole blood in an eppendorf tube, then add 10 μ l of saponin. Saponin is a differential lysing agent. If you lyse blood with only water, both RBC and WBC will be lysed. Saponin perforates only RBC. Vortex the mixture, leave for 3mins. Spin for 3 mins. Decant supernatant.
2. Wash off RBC debris & Hb with 0.05% saponin in normal saline (NaCl) X2cc
3. Lyse WBC by adding
 - Suspend pellet in 25 μ l Neucleo lysing buffer
 - Add 5 μ l Proteinase K (25mg/100 μ l)
4. Incubate at 56 $^{\circ}$ C-2hrs
5. Extract DNA by adding 300 μ l distilled H₂O to the mixture. Vortex. Add equal volume (\approx 500 μ l) of phenol-chloroform mix (Phenol (25) - chloroform (24) - isoamy l (1)). Centrifuge at 13,000 – 15,000 RPM for 5mins. Remove lower layer. Repeat extraction twice.
6. Add 500 μ l of chloroform and vortex. Spin for 5mins and remove lower layer. X2cc
7. Remove the aqueous top layer into a pre-labeled tube and add 2 vols of cold Absolute ethanol (\approx 1ml) mix gently. DNA strand may become visible.
8. Spin and decant supernatant. Wash DNA in 70% EtOH. Spin and decant supernatant.
9. Dissolve DNA in distilled H₂O (DNA can be denatured with time) or TE; which is a chelating agent. it removes ion cofactors like Ca²⁺, Mg²⁺ required for DNAses.

PCR REACTION

Material Supplied

Two lyophilized oligonucleotide primers SC2 and SC5 are supplied. Each tube contains 10 nmoles of the primer. The quantity supplied is sufficient for 400 regular 50ml PCR reaction.

Reconstitution

Stock solution: Add 50ml sterile water to each tube containing the primer.

The 10 nmoles of primer when dissolved in 50ml water will give a solution of 200molar i.e. 200 pmoles/ml.

PCR* reaction (see Appendix for Details)

PCR Reaction mix

	50µl rxn	1ml rxn
10 x PCR Buffer	4.5µl	100µl
dNTP mix (2.5 mM each)	4µl	100µl
Primer Mixz (10 pmal/ml each) (25 pmol of each primer/50m /)	2.5µl	63µl
Sterile water	34µl	737µl
Total	45µl	1ml

β globin primer

SC5

SC2

PCR profile

Denaturation 94°C 30 sec.

Annealing 58°C 30sec.

Elongation 72°C 1 min.

30 cycles, 7min. 72°C extension, 4°C soak.

Restriction enzyme digestion (100m 1 reaction)

DNA 45ml (45ml of 50ml PCR reaction)

10 x buffer 10ml

Dde-I 10-30 units

Sterile water up to 100ml

Precipitate after over night digestion, dissolve pellets in 5ml 1x loading buffer.

A. Electrophoresis

Load samples to 1.5% agarose gel. Run at 90 mAmps for 2.5 hrs.

B. Results

Mutation abolishes restriction site.

Figure 1: PCR gel

PCR product Fragment Size 233 bp

Fragment Sizes After Dde I Digestion

A/A	A/S	S/S
178+55bp	233+178+55bp	233 bp




Figure 1, Typical Sickle cell genotype analysis of PCR product digested with Dde I. Lane 1 is molecular weight markers. Lane 2 is undigested PCR product. Lanes 3, 4 and 6 is DNA with A/S genotype . Lane 5 is A/A genotype DNA and Lane 7 represents DNA with S/S genotype.

References

1. Saiki et al. (1985) Science 230 1350-1354
2. Wu et al. (1989) PNAS 86:2757-2760
3. Conner et al. (1983) PNAS 80:278-282

**The polymerase chain reaction (PCR) process is covered by patents owned by Hoffmarin-La Roche. A license to perform is automatically granted by the use of authorized reagents.

Appendix

PCR Premix preparation

Typical Premix

PCR Reaction mix

	50 μ l rxn	1ml rxn
10 x PCR Buffer	4.5 μ l	100 μ l
dNTP mix (2.5 mM each)	4 μ l	100 μ l
Primer Mixz (10 pmal/ml each) (25 pmol of each primer/50m /)	2.5 μ l	63 μ l
Sterile water	34 μ l	737 μ l
Total	45μl	1ml

Nucleotide Dilution

Stock-----100 mM: Prepare a final diluted 2.5 mM solution

Preparation

Each 100 mM dNTP	25 μ l (Total 500 μ l)
Water	4.5 ml

Total Volume 5.0ml

Taq Premix (per 50 μ l reaction, scale up as required)

10 x PCR Buffer	0.5 μ l
Taq polymerase	0.25 μ l
Sterile water	4.25 μ l
Total	5μl/rxn

PCR reaction (50 μ l)

Diluted DNA (100ng/ml)	1 μ l
PCR premix	4.5 μ l
Taq premix	5 μ l

PCR products post-processing

1. For all layered PCR only Acid 200ml of CHC13 to each tube, vortex and spin

2. Transfer the upper aqueous layer to a fresh eppendorf tube, add 1/10 volume of 3M NaAc (pH 5.2) and 2 volumes of absolute ethanol, precipitate DNA at -80°C for 10 minutes.
3. Spin, rinse the DNA pellet with 700ml of 75% ethanol and dry the pellet in the speedvac
4. Dissolve the pellet in adequate amount of TE.